

User Guide


Voronoi Health Analytics Inc.







Home

This is an interactive user guide. Click on the blue boxes to learn more about each part of DAFS or click on a menu item on the right.


You can return here anytime to check out other features by clicking the  icon.



Use these buttons to open helpful menus!



Overview






Preprocessing




Preprocess

Preprocessing will detect all DICOMs within an Input Folder and use the embedded DICOM metadata information to *sort* and *copy* each unique scan into its own folder. It is integral to perform Preprocessing before proceeding with any other processing.



Input Folder – Select any top-level folder containing DICOMs.

- All DICOMs, including those in subfolders within the designated Input Folder will be preprocessed.
- Ensure that no DAFS System Folders ([what's this?](#)) are in the Input Folder path.



Output Folder – Select any folder. DAFS will create a new folder known as the Working Folder ([what's this?](#)) within the Output Folder (date & time stamped) and the curated scans will be inside, copied from the Input Folder.

- DAFS will create a Working Folder inside the Output Folder.
- Ensure that no DAFS System Folders are in the Output Folder path.



Disable strict scan integrity checking – DAFS has strict checking for the integrity of the 3D scan reconstructed from the DICOMs. Enable this to override strict checking for slice orthogonality, heterogeneous acquisition numbers and heterogeneous slice orientations.

Preprocessing may take some time depending on the size of the Input folder. For an Input folder containing ~10000 files, curation will take ~12.5 minutes.

Where can I find out which scans did not get preprocessed?

[Learn more](#)



Preprocess: Configuring CT/PT Mapping Rules

CT and PET mapping rules can be configured in:

%LOCALAPPDATA%\Voronoi Health Analytics\DAFS\dafs_config.json

Example:

```
{  
  "version": "0.0.2",  
  "match_attrs": ["PatientID", "StudyInstanceUID"],  
  "match_regex": {"SeriesDescription": ".*[^N]AC.*"}  
}
```

"match_attrs": Mapping based on matching attributes: Specify a list of attributes that must be exactly matching between the CT and PET.

"match_regex": Regex matching: Specify a list of attributes which values must match specific regular expressions. We recommend using this in conjunction with "match_attrs", not on its own.

After Preprocessing is run and there are both valid CT scans and PT scans in your dataset, a CT to PT mapping csv will be made in the Working Folder with matching CT-PT pairs according to the mapping rules. This csv can be manually edited if required.



Processing and Exporting



Working Folder

A Working Folder is generated within the Output Folder following Preprocessing and is the location in which DAFS and CAST will both read and write. A Working Folder will be generated with a date & time stamp. DAFS will update the Working Folder each time Start Processing is clicked.



Working Folder – Select a folder generated by Preprocessing. This folder will likely have a date and time stamp if the name has not been modified. Generated after preprocessing is complete.

- DAFS and CAST will read and write to this folder.

What's inside a Working Folder?

[Learn more](#)

Click to see a full map of a Working Folder.

Curated Scans

These are folders, each containing a unique scan that has been curated via Preprocessing.

DAFS-report_date_time folders

These are folders that store processing information and reports from each DAFS run. A DAFS-report folder that is date and time stamped will be generated each time DAFS is run on the same Working Folder, acting as a historical record of DAFS processes. Each DAFS-report may include one or more of the following subfolders depending on the options chosen.



logs

Contains logs on errors that may have occurred during processing. The contents in this folder may help with troubleshooting if the DAFS output was not as expected. Generated every run - **It is highly recommended to check this folder after processing.**



CSVs

Contains results from Preprocessing, Anonymization and Publish in csv format.



quickchecks

Contains visual media to aid with quality assessment in PNG format.



reports

Contains results from enabling PDF reports



Representation of a folder structure after three DAFS runs following a suggested workflow.

For more suggested workflows...

[Learn more](#)

Output Folder

Working Folder

Preprocessing must always be run first.

1st run DAFS-report
• Preprocessing

csvs
DICOM Curation.csv
logs
*****-failures-log.txt*

Find the curated ID to original file path map here.

2nd run DAFS-report
• Anonymize
• Segment Tissues
• Annotate slices
• Quickchecks

DAFS-report_2022-00-01_00.00.00.000
logs
******-failures-log.txt*
quickchecks
csvs
Anonymization csv
reports

Any processing failures can be found in the logs folder.

Find quickchecks here for quality assessment.

Find an anonymization report here.

You may also choose to directly Publish without quality control in this run.

Find PDF reports here.

3rd run DAFS-report
• Publish

DAFS-report_2022-00-02_00.00.00.000
logs
******-failures-log.txt*
csvs
Publish csvs

Find measurement results here for all scans in the Working Folder.

You may choose to run Publish and Quickchecks in multiple runs to generate the most up-to-date results after segmentation and/or annotation quality control.

scan_1
scan_2
scan_3
scan_4
scan_5
scan_6

These are folders each containing a unique, curated scan generated by Preprocessing in the 1st run.



Auto-Label

Enable automated segmentation and landmarking.

Reset Labels

Enabling this checkbox will prompt DAFS to replace existing segments. If manual corrections have been done in CAST, it is likely that you should NOT use this feature. Segmentations that are replaced cannot be recovered. **Please use with caution.**



Quickchecks

Generate visual media for efficient quality checking of segmentations in the Working Folder. Run before, during and/or after segmentation quality control to update quickchecks with the most recent changes. Find results in the quickchecks folder.

Entries

Slices

Head;C1;C2;C3;C4;C5;C6;C7;T1;T2;T3;T4;T5;T6;T7;T8;T9;T10;T11;T12;L1;L2;L3;L4;L5;L6;LSTV;Sacrum;Femur;Knee;Tibia;Feet

Make Quickchecks with slice levels selected in the text area. Defining Slices will create sagittal and coronal .pngs that display the annotated slices in the image. Axial images will also be sampled from selected Slice levels. If Slices are selected in conjunction with Tissues, Slices will take priority in axial quickcheck sampling. Leave empty for equal sampling throughout the scan.

Tissues [Tissue Dictionary](#)

Make Quickchecks with Tissues selected in the text area. Using tissues will create sagittal and coronal .pngs that are sampled within the ROI of the tissue. Axial images will also be sampled within the ROI of the tissue **if no Slices are selected**. If Slices are selected in conjunction with Tissues, Slices will take priority in axial quickcheck sampling.

Opacity Control

Control the level of opacity level of segmentation displayed on quickchecks.

Identifier

Enable printing of an image identifier on the quickcheck.

- Options include None, PatientID, PatientName, SeriesNumber, StudyID, AccessionNumber and MRN.
- If the selected identifier does not exist for a scan, no identifier will be printed on the quickcheck.

PNG files naming convention

- annotation-coronal/annotation-sagittal – images displaying annotated slices.
- tissues_and_annotations-axial – images displaying tissues and annotated slices, with annotated slices taking priority in sampling if selected.
- tissues-coronal/tissues-sagittal – images displaying tissues.



Anonymize

Find the anonymization report in the csvs folder.

<u>Attribute</u>	<u>Method of anonymization</u>
InstanceCreatorUID	remove
InstitutionAddress	remove
ReferringPhysicianAddress	remove
ReferringPhysicianTelephoneNumbers	remove
StationName	remove
StudyDescription	remove
SeriesDescription	remove
InstitutionalDepartmentName	remove
PhysiciansOfRecord	remove
NameOfPhysiciansReadingStudy	remove
OperatorsName	remove
AdmittingDiagnosesDescription	remove
DerivationDescription	remove
PatientBirthTime	remove
OtherPatientIDs	remove
OtherPatientNames	remove
PatientAge	remove
PatientSize	remove
PatientWeight	remove
PatientAddress	remove
MedicalRecordLocator	remove
EthnicGroup	remove
Occupation	remove
AdditionalPatientHistory	remove
PatientComments	remove
ProtocolName	remove
ImageComments	remove
RequestAttributesSequence	remove
ContentSequence	remove
StorageMediaFileSetUID	remove

<u>Attribute</u>	<u>Method of anonymization</u>
SOPInstanceUID	change
AccessionNumber	change
InstitutionName	change
ReferringPhysicianName	change
PerformingPhysicianName	change
ReferencedSOPInstanceUID	change
PatientName	change
PatientID	change
PatientBirthDate	change
PatientSex	change
StudyInstanceUID	change
StudyID	change
FrameOfReferenceUID	change
SynchronizationFrameOfReferenceUID	change
ID	change
UID	change
ReferencedFrameOfReferenceUID	change
RelatedFrameOfReferenceUID	change

Non-anonymizable attributes (3D scan integrity)

SeriesInstanceUID	PixelSpacing
AcquisitionNumber	SliceLocation
SeriesNumber	SliceThickness
AnatomicalOrientationType	RescaleIntercept
ImageOrientationPatient	RescaleSlope
ImagePositionPatient	PhotometricInterpretation
PixelData	BitsAllocated
Rows	BitsStored
Columns	HighBit



Publish

Compile information derived from DICOMs, segmentations and landmarks into csv files. Multiple entries in textboxes are possible by delimiting each entry with a semicolon. See the [Publish CSV Dictionary](#) for column definitions. **Cross sectional areas for volumetric slabs are averaged across all slices within the slab.*

Region-Of-Interest (ROI) Constructor

Derive from Slices and Slab(s)

Use anatomical locations to derive tissue measurements by defining vertebral level Slice Annotations (inputs are case insensitive).

Examples:

- **FULL_SCAN** - use measurements from the entire scan.
- **L3mid** - use measurements from L3mid slice only (cross-sectional area only).
- **T1start-to-L5end** – use measurements from the top of T1 vertebrae to bottom of L5 vertebrae.
- **avg-T12mid[3]** – use measurements at T12mid averaged across 3 slices above and below.
- **L1start-to-L1end** – use measurements from the top of L1 to bottom of L1 vertebrae.
- **L1** – equivalent to L1start-to-L1end.
- **ScanStart** – defined as the top of the scan; the first slice. ScanEnd is defined as the bottom of the scan; the last slice.

Derive from Tissues and Parts of Tissues

[Tissue Dictionary](#)

Select tissues and HU ranges to derive measurements from.

Examples:

- **ALLSKM** – derive measurements from total skeletal muscle within the scan.
- **ALLSKM[-29,150]** – derive measurements from total skeletal muscle within the range of -29 to 150 HU only found within the scan.
- **VAT[-150,-50]-U-SAT[-190,-30]** – derive measurements from VAT within -150 to -50 HU and SAT within -190 to -30 HU combined (unionized) as one segment. Multiple ROIs can be unionized separated by -U-.
- **LIV** – derive measurements from liver segmentation.
- **SPL** – derive measurements from spleen segmentation.
- **LKID** – derive measurements from the left kidney segmentation.

Print DICOM Attributes – Choose DICOM header information to include in the CSV

Define a set of DICOM attributes to be printed as columns in the csv. This is useful if you wish to sort scans by certain DICOM attributes when viewing the csv report. Anonymized attributes will print as NA.

How do I view EUR format in MS Excel?

[Learn more](#)

My data has leading zeroes and is not being displayed correctly in Excel.

[Learn more](#)



Tissue Dictionary p1

Muscle Group

Colour	Entry	Anatomical ROI
	ALLSKM	Total skeletal muscle within the scan
	ALLSKM_NOARMS	Total skeletal muscle excluding arms
	LPECMJR	Left pectoralis major
	RPECMJR	Right pectoralis major
	LPECMNR	Left pectoralis minor
	RPECMNR	Right-pectoralis-minor
	LTEMPORALIS	Left temporalis
	RTEMPORALIS	Right temporalis
	LMASSETER	Left masseter
	RMASSETER	Right masseter
	LILIOPSOAS	Left iliopsoas
	RILIOPSOAS	Right iliopsoas
	LUPLGSKM	Left upper leg skeletal muscle
	RUPLGSKM	Right upper leg skeletal muscle
	LLWLGSKM	Left lower leg skeletal muscle
	RLWLGSKM	Right lower leg skeletal muscle
	LASKM	Left arm skeletal muscle
	RASKM	Right arm skeletal muscle

Ventral Adipose Tissue Group

Colour	Entry	Anatomical ROI
	VAT	Visceral adipose tissue
	EpAT	Epicardial adipose tissue
	PaAT	Paracardial adipose tissue
	ThAT	Thoracic adipose tissue

IMAT (Intra-muscular Adipose Tissue) Group

Colour	Entry	Anatomical ROI
	ALLIMAT	Intramuscular adipose tissue
	ALLIMAT_NOARMS	Total scan IMAT excluding the arms
	LUPLGIMAT	Left upper leg IMAT
	RUPLGIMAT	Right upper leg IMAT
	LLWLGIMAT	Left lower leg IMAT
	RLWLGIMAT	Right lower leg IMAT
	LAIMAT	Left arm IMAT
	RAIMAT	Right arm IMAT

SAT (Subcutaneous Adipose Tissue) Group

Colour	Entry	Anatomical ROI
	SAT	Subcutaneous adipose tissue
	LASAT	Left arm SAT
	RASAT	Right arm SAT

More tissues on the next page



Tissue Dictionary p2

Fluid Buildup		
Colour	Entry	Anatomical ROI
	ASCITES	Ascites; abdominal fluid
	LUNGEFFUSION	Lung effusion; pleural effusion

Cardiovascular		
Colour	Entry	Anatomical ROI
	AOC-U-CAAC	enables AOC-CAAC-based agatston score
	AOC	enables AOC-based agatston score
	CAAC	enables CAAC-based agatston score

Organs					
Colour	Entry	Anatomical ROI	Colour	Entry	Anatomical ROI
	LIV	Liver		AO	Aorta
	SPL	Spleen		AOC	Aortic calcification
	LKID	Left kidney		CAAC	Cardiac aggregate calcium
	RKID	Right kidney		LLUNG	Left lung
	GB	Gallbladder		RLUNG	Right lung
	HRT	Heart		TRACH	Trachea
	PANC	Pancreas		BLADDER	Bladder
	BRAIN	Brain		PROSTATE	Prostate

Bone Group		
Colour	Entry	Anatomical ROI
	ALLBONE	Total bone within the scan
	ALLBONE_NOARMS	Total bone within the scan excluding the arms
	LFMRHEAD	Left femoral head
	RFMRHEAD	Right femoral head
	LFMRNECK	Left femoral neck
	RFMRNECK	Right femoral neck
	LFMRSHAFT	Left femoral shaft
	RFMRSHAFT	Right femoral shaft
	LHPBONE	Left hip bone
	RHPBONE	Right hip bone
	TRBCLR	Trabecular bone
	LABONE	Left arm bone
	RABONE	Right arm bone



Reports

Body Composition Reports

Ensure that you have your annotations ready for your scans before running PDF report generation.

Outputs will be sent to a DAFS Report folder in your working folder under the folder name “reports”.

Specify one or more landmarks with each entry separated by a semicolon. The entries are case insensitive.

Examples:

- `L3mid` (the middle slice of L3).
- `avg-L3mid[5]` (average cross sectional area 5 slices above and below L3 middle annotation).
- `T4start` (the top slice marked as T4).
- `L5end` (the bottom slice marked as L5).
- `T12mid;L3mid;avg-L5mid[2]` (print reports for T12mid, L3mid and L5mid averaged across 2 slices above and below).
- `FULL_SCAN` - use measurements from the entire scan.
- `T1start-to-L5end` – use measurements from the top of T1 vertebrae to bottom of L5 vertebrae.
- `L1start-to-L1end` – use measurements from the top of L1 to bottom of L1 vertebrae.
- `L1` – equivalent to L1start-to-L1end.
- `L4start-to-scanend` – use measurements from the top of L4 to bottom of scan. “scanend” is used to define the lowest slice of the scan.
- `L4start-to-scanstart` – use measurements from the top of L4 to top of scan. “scanstart” is used to define the highest slice of the scan.

Organ Reports

Outputs will be sent to a DAFS Report folder in your working folder under the folder name “reports”.

Specify one or more organs with each entry separated by a semicolon. The entries are case insensitive.

Possible entries:

`LIV;SPL;LKID;RKID;HRT;CAAC;GB;AO;AOC;LLUNG;RLUNG;TRACH;PANC;BRAIN;BLADDER;PROSTATE`



CAST: Segmenter

Select the Working Directory

Toggle blending of the CT/PT (if both are available)
100% CT | 50/50 CT/PT | 100% PT

Toggle crosshair (hover on viewport while pressing shift)

Save changes

Toggle maximized viewport

Link dual views to scroll and zoom together

CT Windows:

Blend Volumes: CT PT

Viewport: Triple Side-by-Side

Center the Volume

Toggle slice intersections

Toggle segmentation outline

Change editor

Change windows of the CT volume

Choose the layout of the Viewport

Center the Volume

Toggle slice intersections

Toggle segmentation outline

Change editor

List of scans in the Working Directory; Double-click to load

- scan_1_2024-04-25_11.29.54.067
- scan_2_2024-04-25_11.29.54.067
- scan_3_2024-04-25_11.58.41.237
- scan_4_2024-04-25_11.58.41.237

Load

Scan Information

Scan: scan_3_2024-04-25_11.58.41.237 Patient ID: ANON41119

Patient Name: DOE JANE Study ID: 2.16.840.1.114362.1.12177026.2...

Patient Age: 041Y Series ID:

Patient Sex: Seri

KVM: Vox

MA: Set the editable intensity ranges

Colour	Tissue	Lock	Hide	Min HU	Max HU	Reset HU
SKM	SKM	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	-32768	32767	<input type="button" value="↺"/>
RPECMJR	RPECMJR	<input type="checkbox"/>	<input checked="" type="checkbox"/>	-32768	32767	<input type="button" value="↺"/>
LPECMJR	LPECMJR	<input type="checkbox"/>	<input checked="" type="checkbox"/>	-32768	32767	<input type="button" value="↺"/>
RPECMNR	RPECMNR	<input type="checkbox"/>	<input checked="" type="checkbox"/>	-32768	32767	<input type="button" value="↺"/>
LPECMNR	LPECMNR	<input type="checkbox"/>	<input checked="" type="checkbox"/>	-32768	32767	<input type="button" value="↺"/>
LTEMPORA...	LTEMPORA...	<input type="checkbox"/>	<input checked="" type="checkbox"/>	-32768	32767	<input type="button" value="↺"/>
RTEMPORA...	RTEMPORA...	<input type="checkbox"/>	<input checked="" type="checkbox"/>	-32768	32767	<input type="button" value="↺"/>

Filter: Unfiltered

Fill Opacity:

Equip: Tool: None 3D Brush Off

Diameter:

DICOM Information panel for the loaded scan

List of segmentations; click to enable editing.

Filter for segments in the Segment list. Use "ALL" to show all segments

Adjust segmentation fill opacity

Change editing tool; Shortcuts: 1, 2, 3, 4 Pen, Erase, Erase All, Pointer

Lock and/or hide filtered out segments (prevent overwriting and visualization)

Toggle 3D brush

Adjust brush diameter









Filtering the Segment list

Type in tissue names separated by semi-colon or hyphen and press Enter to display them in the Segment list, thereby, narrowing the focus of segmentation QC to only segments you wish to see and edit. Examples:

















- ALL will display all available segments for the scan.
- SKM;VAT will display SKM and VAT segments only.
- SKM-VAT will display all segments between SKM and VAT.

For Unfiltered segments:

-   Segments not in the Segment list cannot be overwritten or seen.
-   Segments not in the Segment list can be overwritten but not seen.
-   Segments not in the Segment list can be overwritten and seen.

Keyboard and mouse shortcuts

(* requires mouse cursor to be on Viewport)

			Pen	Erase	Erase All	Pointer
Toggle panel on/off	 + 	Select tool				
Zoom*	 or  +  scroll	Cross reference	 + 			
Pan*	 or  + 	Hide/show segmentation				
Apply or erase (tool dependent)*						



CAST: Annotator

Annotator is a unique tool that enables editing of annotations generated by AVA or creation of manual annotations without AI assistance. While there is a small learning curve, the next three pages will help you get familiar with Annotator and its functions.

DICOM Information panel for the loaded scan


The screenshot shows the 'Scan Information' panel with a 'Mark Slab Endpoints' section. It features a grid of 28 buttons for anatomical landmarks: Head, C1-C7, T1-T12, L1-L6, Sacrum, LSTV, Femur, Knee, Tibia, and Feet. Above the grid are 'Mark as Middle Slice' and 'Unmark Middle Slice' buttons. Below the grid are 'Delete all', 'Delete', 'Redo', and 'Undo' buttons.

Mark or unmark a slice as the "middle" of the landmark

Buttons to add or modify slab endpoints (heights of slabs)

Deletes all annotations

Deletes the slab that corresponds with the current axial slice shown in the viewport




Annotator enables assignment and editing of vertebral annotations to DICOM slices. The tool works by adjusting landmark “slab” heights. Slab heights can range from 1 slice to N number of slices in the image. The -start and -end indicate the top and bottom boundaries of the slab and every slice between the two boundaries are annotated as the same landmark. There is also a -mid suffix which indicates the slice annotated as a middle slice of that landmark; this annotation is displayed in white in the ViewPort.

To begin, first try placing your mouse cursor on the sagittal viewport. Hold SHIFT + move your mouse up and down. Notice that the axial view intersection line is following your cursor and that the axial view is changing rapidly.


Next, place the axial view intersection line a few slices above a –start boundary. Then, click the button of the landmark you are working on. Notice that the –start boundary has moved and the height of the slab is increased. Similarly, placing the intersection line below the boundary will decrease the height and contract the slab. Try this with the –end boundary, but this time placing above the boundary will contract slab, and below will expand the slab.

Annotator will always edit the boundary that is closest to the Axial view intersection line.





To assign middle slices, both boundaries for a given landmark must be present. You cannot assign a middle slice on a slab with thickness = 1 (only 1 DICOM is annotated). To reassign a different DICOM as the middle slice, simply place axial slice intersection line within a landmark's slab boundary and click "Mark as middle slice".



You may also mark –start and –end slices as middle slices. Notice that the boundary marked as middle will appear slightly brighter.



You may choose to unmark a middle slice for a given vertebra. Locating the axial view intersection line within the landmark's slab will enable you to reassign slices to the middle slice or remove the middle slice annotation for that vertebrae.



Appendix



Preprocessing Error Logs in Detail

Preprocessing failures will occur if the following errors are detected. You can find a curation failures log in the DAFS-report and logs folder. See [page 6](#) to learn more about where to find reports and logs.


- heterogeneous_series_numbers – Corrupted or empty DICOM header.
- heterogeneous_acquisition_numbers – Corrupted or empty DICOM header.
- heterogeneous_series_uids – Corrupted or empty DICOM header.
- invalid_anatomy_type – Corrupted or empty DICOM header.
- invalid_image_orientation_patient – Corrupted or empty DICOM header.
- invalid_image_position_patient – Corrupted or empty DICOM header.
- heterogeneous_slice_orientations – Corrupted or empty DICOM header.
- duplicate_slice – Corrupted scan containing 2 or more of the same DICOM file.
- no_valid_dicom_slices - Only .asg files are present in the Input folder.
- non_anonymized – PHI fields detected.
- non_axial_orthogonal_orientation – There are non-axial slices present.
- missing_slices – There are DICOMs missing and disrupting 3D scan integrity.
- 3D_scan_loading_failure – Issue with 3D scan.

[Go back to Preprocessing](#)



Generating EUR csv format


Windows Control Panel -> Clock and Region -> Region



Open Excel, and click Open -> Browse -> select All files (*.*) and navigate to the Output folder 'csvs' subfolder to find the *Publish.txt* and open.



A Text import Wizard will appear.



1. Choose file type that best describes your data. It should say Delimited by default. Click Next.


2. Select Comma as the only Delimiter. Click Next.

3. Click Finish.

Go back to Publish




Suggested Workflow



[Go back to Working Folder](#)



Folders Dictionary



Input Folder – Select any top-level folder containing DICOMs. All DICOMs, including those in subfolders within the designated Input Folder will be processed. Applies to Preprocessing only.



Output Folder – Select any folder. DAFS will create a Working Folder within the Output Folder. Applies to Preprocessing only.



Working Folder – Select any curated folder generated by Preprocessing. This folder will likely have a date and time stamp if the name has not been modified and must only have subfolders inside. All processing by DAFS and edits made in CAST will be saved within the Working Folder. DAFS will automatically designate the Working Folder after Preprocessing is complete.



DAFS System Folders – Folder names that cannot be anywhere in the path of the Input Folder, Output Folder or Working Folder. DAFS will not continue processing if the following keywords are detected in any of these paths (case insensitive). Furthermore, any of these key words followed by an underscore (_) cannot be anywhere in the paths.

- DAFS-report
- logs
- quickchecks
- reports
- csvs
- tmp
- Anonymize
- Tissue Statistics_slicewise
- Tissue Statistics
- Preprocessing
- ABACS 3D
- AVA

[Go back to Preprocessing](#)

[Go back to Working Folder](#)



Supported GPUs

GPU utilization can be enabled for automated segmentation and annotation processing. Using GPUs may decrease processing time by 10 times or more when compared to CPU.

Recommended GPU: Nvidia GeForce RTX 2070 or higher.

List of supported Nvidia GPU microarchitectures and tested graphics card models:

- Ada Lovelace
 - Nvidia GeForce 4000 series
- Ampere
 - Nvidia GeForce RTX 3060
 - Nvidia GeForce RTX 3060 Ti
 - Nvidia GeForce RTX 3070
 - Nvidia GeForce RTX 3090
 - Nvidia RTX A5000
- Turing
 - Nvidia GeForce RTX 2060
 - Nvidia GeForce RTX 2070
 - Nvidia GeForce RTX 2070 Super
- Volta
- Pascal
 - Nvidia GeForce MX150
 - Nvidia GeForce MX250
 - Nvidia GeForce MX330
 - Nvidia GeForce GTX 1070
 - Nvidia GeForce GTX 1050 Ti

Note: DAFS 3 may also be compatible with other graphics card models not listed above with the supported microarchitectures.

If you think you may have a compatible GPU, but DAFS is not detecting it, please contact us at support@voronoihealthanalytics.com.



Publish CSV Dictionary

Category General format

General scan attributes: `scan_folder, orig_path, voxel_size_WxLxH`

Slab attributes: `$slab;num_slices, $slab;img_size_WxLxH, $slab;height_mm`

Slab and Tissue (ROI) attributes: `$slab;$tissue;$tissue_statistic`

- `$tissue_statistic` = `cross.sectional_area_pixels, cross_sectional_area_cm2, volume_voxels, volume_cm3, HU_mean, HU_std, HU_min, HU_max`

** "\$" represents a variable.**

Term	Meaning
<code>defined_slab</code>	(applicable to <code>by_slab</code> csv only) Describes the slab in which the measurements belong.
<code>scan_folder</code>	A unique ID assigned to a scan after curation which also is the folder name of that scan.
<code>orig_path</code>	The original path of the scan or DICOM files.
<code>img_size_WxLxH</code>	Image dimensions in width, length and height.
<code>voxel_size_WxLxH</code>	Voxel dimensions in width, length and height.
<code>height_mm</code>	The vertical anatomical height of the defined.slabs in mm.
<code>num_slices</code>	The number of DICOM slices in the defined.slabs.
<code>cross.sectional_area_pixels</code>	Cross sectional number of pixels - number of voxels within the region-of-interest.
<code>cross_sectional_area_cm2</code>	Cross sectional area in centimeters squared within the region-of-interest.
<code>volume_voxels</code>	Volume number of voxels - number of voxels within the region-of-interest.
<code>volume_cm3</code>	Volume in centimeters cubed within the region-of-interest.
<code>HU_mean</code>	Mean Hounsfield Unit within the region-of-interest.
<code>HU_std</code>	Hounsfield Unit standard deviation within the region-of-interest.
<code>HU_min</code>	Minimum Hounsfield Unit within the region-of-interest.
<code>HU_max</code>	Maximum Hounsfield Unit within the region-of-interest.
<code>pet_uptake_*</code>	Measurements derived from the PET image
<code>tissue[#####]</code>	The region-of-interest defined by the tissue name and hounsfield unit range.
<code>tissue[#####]-U-tissue[#####]</code>	An ROI defined by multiple region-of-interests unionized defined by the tissue names and Hounsfield unit ranges.

[Go back to Publish](#)



License and Information

Expiration date – Date of license expiry.

Scans processed – The number of unique scans that have been processed out of the total number of available scans provisioned in the license.

Scans remaining – The number of available scans left in the license.